

## EXPERT PLUS / MULTIREAD 400 PRINTERS

### 2011 COMPATIBLE MODELS



The following list of printers are specified as being directly compatible or with the use of language emulation by their respective manufacturers.

Please note that most printers require additional configuration in order to accept the printer languages output from our instruments. This can sometimes be set via a connected PC to the printer or from the printer keypad and display (if fitted). The printer should be set in emulation mode to accept one of the following languages: PCL4, PCL5, Epson FX, Epson ESC/P or ESC/P2. Please consult the printer manufacturer's user manual for correct setup instructions.

The following printers should be selected from the printer selection menu on our instruments screen depending on the emulated language set up on your printer:

Printer	Emulated Language on set on printer
Compatible > LaserJet 4	PCL4, PCL5
Epson > FX-850	Epson FX
Compatible > Epson ESC/P	Epson ESC/P
Compatible > Epson ESC/P2	Epson ESC/P2

Please contact us if you require further assistance.

#### Compatible Printers - 2011

- Dell 2350d
- Dell 2350dn
- Epson M2000D
- Epson M2000DN
- Epson M2000DTN
- Epson M2000DT
- Epson M2300D
- Epson M2300DN
- Epson M2300DT
- Epson M2300DTN
- Epson M2400D
- Epson M2400DT
- Epson M2400DTN
- Epson EPL-6200
- Epson EPL-6200N
- HP LaserJet P2035
- Kyocera FS-2020D
- Lexmark E260
- Lexmark E260d
- Lexmark E260dn
- Lexmark E360
- Lexmark E360d
- Lexmark E360dn
- OKI B430DN
- OKI B410d
- OKI B411
- OKI B411d
- OKI B411dn
- OKI B431
- OKI B431dn
- OKI B440dn
- OKI B2400
- OKI B2400n
- OKI B4600
- OKI B4600n
- OKI B4600PS
- OKI B4600nPS
- Samsung ML-3471ND
- Xerox Phaser 3600