


# UVDAS Accessories - Low volume cells



80-2103-69



80-3000-81



80-2103-68



80-3000-83

- The Biowavell/Lightwavell/Biowave DNA range are compatible with
  - 70ul Quartz cell(80-2103-69)
  - 70ul disposable UV Plastic cell (80-3000-81)
  - 10ul 10mm pathlength Quartz cell (80-3000-83)
  - 7ul 5mm pathlength Quartz cell (80-2103-68) note - external dimensions are same as standard 10mm pathlength cell